


# ORL


ORL drape-universal, reinforced, with adhesive aperture.

<b>Code</b>	106-102-S
<b>Material</b>	2SBL
<b>Size (cm)</b>	200 x 240
<b>Reinforcement (cm)</b>	Vis28 60 x 75
<b>Aperture (cm)</b>	7
<b>Accessories</b>	AT



ORL drape, with adhesive aperture.

<b>Code</b>	106-103-S
<b>Material</b>	DS2
<b>Size (cm)</b>	200 x 240
<b>Reinforcement (cm)</b>	-
<b>Aperture (cm)</b>	7
<b>Accessories</b>	AT



Tracheotomy drape reinforced, with adhesive aperture.

<b>Code</b>	501-008-S
<b>Material</b>	SMS
<b>Size (cm)</b>	160 x 300
<b>Reinforcement (cm)</b>	DS2 75 x 100
<b>Aperture (cm)</b>	10
<b>Accessories</b>	AT

